

Chadda Surgicals - Product List

Cell: Kawal Chadda : +919819691699, Mohit Chadda : +919930238882

Email: chaddasurgicals@hotmail.com

LATEST TITANIUM MINI BONE PLATES & SCREWS

- ❖ BASIC INSTRUMENT SET FOR MINI PLATES IN S.S. BOX (THE SET CONSISTS OF PLATE BENDERS, SCREW DRIVERS, FORCEPS FOR HOLDING BONE PLATE AND SCREW, SCALE, GAUGE) TITANIUM
- ❖ IMPLANTS SET OF MINISCREWS AND MINI PLATES IN S.S. BOX
 - 1.5MM, 2MM OR 2.5MM DIA ASSORTED SCREWS (6 - 12 MM) 100 NOS. TITANIUM
 - MINI PLATES - ASSORTED (UP TO 6 HOLES) 25 NOS. TITANIUM
- ❖ MINI SCREW- 1.5MM DIA SLOT HEAD (2MM) LENGTH 4,6,8,10 MM.
- ❖ MINI SCREW- 2MM DIA SLOT HEAD (3MM) LENGTH 4,6,8,10,12,14,16 MM.
- ❖ MINI SCREW- 2.5MM DIA SLOT HEAD (4MM) LENGTH 4,6,8,10,12,14,16 MM.
- ❖ MINI PLATES, TITANIUM - 1.5, 2, 2.5 MM OR SCREW SIZE :
 - 2 HOLES- CONTINUOUS OR WITH GAP
 - 3 HOLES- CONTINUOUS OR WITH GAP
 - 4 HOLES- CONTINUOUS OR WITH GAP
 - 5 HOLES- CONTINUOUS OR WITH GAP
 - HOLES- CONTINUOUS OR WITH GAP
 - 16 HOLES- CONTINUOUS OR WITH GAPFOR EVERY ADDITIONAL HOLE
 - ORBITA- 4 HOLES.
 - ORBITA- 6 HOLES.
 - 'L' SHAPED- SMALL- 4 HOLES RIGHT OR LEFT.
 - 'L' SHAPED- BIG- WITH GAP - 4 HOLES RIGHT OR LEFT, BIG WITH GAP3 X 3 RIGHT OR LEFT
 - 'T' SHAPED- BIG 3 X 3 HOLES
 - 'L' SHAPED 6 HOLES
 - 'Y' SHAPED 6 HOLES
 - 'X' SHAPED 6 HOLES
 - 3D PLATE- 4 HOLES
- ❖ RECONSTRUCTION PLATES FOR MANDIBLE FOR 1.7 OR 2.5 MM SCREWS
 - UPTO 20 HOLES, TITANIUM
 - STRAIGHT
 - ANGLED RIGHT OR LEFT
 - FOR EVERY ADDITIONAL HOLE
- ❖ LAG SCREW, S.S. :
 - 2 MM OR 2.5 MM DIA
 - 8 MM TO 20 MM LONG
 - 2.5 MM DIA - 25, 30 MM LONG
- ❖ EMERGENCY SCREW, S.S. :
2, 2.2, 2.5, 2.7 MM THREAD, DIA UP TO 25 MM LONG.
- ❖ TMJ CONDYLE PLATES, TITANIUM : RIGHT, LEFT. SMALL OR LARGE.
- ❖ TMJ CONDYLE WITH RECONSTRUCTION PLATE, TITANIUM.

Chadda Surgicals - Product List

Cell: Kawal Chadda : +919819691699, Mohit Chadda : +919930238882

Email: chaddasurgicals@hotmail.com

- ❖ K WIRES, S.S. 0.5, 1, 1.5 MM- 3", 4" LONG.
- ❖ S.S. WIRE REEL - 24, 26, 28, 30 GAUGE.
- ❖ ARCH BAR 12" LONG, S.S.
- ❖ DRILL BITS, TITANIUM - 1, 1.5, 2 MM DIA.
- ❖ HEISTER JAW OPENER, S.S.
- ❖ SELF HOLDING SCREW DRIVER, 1.5 MM, 2 MM, 2.5 MM, S.S.
- ❖ SELF HOLDING SCREW DRIVER, 1.5 MM, 2 MM, 2.5 MM, TITANIUM.
- ❖ PLAIN SCREW DRIVER, 1.5 MM, 2 MM, 2.5 MM, S.S.
- ❖ PLAIN SCREW DRIVER, 1.5 MM, 2 MM, 2.5 MM, TITANIUM.
- ❖ ELEVATOR STRAIGHT, S.S.
- ❖ COPLAND ELEVATOR, S.S.
- ❖ PLATE BENDING PLIERS WITH PIN, S.S.
- ❖ MOULDING LEVER, S.S.
- ❖ BONE HOLDING CLAMP, S.S.
- ❖ DEPTH GAUGE, S.S.
- ❖ SCREW/PLATE HOLDING FORCEPS, S.S.
- ❖ VERTICAL FIXATOR, S.S.
- ❖ HAND DRILL, SMALL, MANUAL, S.S.